

Product Datasheet

Human PLGF-2 Protein (orb1475249)

Description	Recombinant Human Placenta Growth Factor is produced by our Mammalian expression system and the targ...
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB,150mM NaCl ,PH7.2.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	P49763
Expiration Date	6 months from date of receipt.