

## Product Datasheet

### Human IL-17RB Protein (orb1475315)

<b>Description</b>	Recombinant Human Interleukin-17 Receptor B is produced by our Mammalian expression system and the t...
<b>Endotoxins</b>	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9NRM6</b>
<b>Expiration Date</b>	6 months from date of receipt.