

Product Datasheet

Human Siglec-6 Protein (orb1475320)

Description	Recombinant Human Sialic acid-binding Ig-like lectin 6 is produced by our Mammalian expression syste...
Endotoxins	Less than 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 μ m filtered solution of 20mM Tris-HCl,150mM NaCl,pH 8.5.
Storage	Lyophilized protein should be stored at \leq -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at \leq -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Uniprot ID	O43699
Expiration Date	6 months from date of receipt.