

Product Datasheet

Human MAdCAM-1 Protein (orb1475321)

Description	Recombinant Human MADCAM1 is produced by our Mammalian expression system and the target gene encodin...
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	Q4PKD0
Expiration Date	6 months from date of receipt.