

Product Datasheet

Mouse PVRIg Protein (orb1475327)

| | |
|------------------------|---|
| Description | Recombinant Mouse Transmembrane Protein PVRIg is produced by our Mammalian expression system and the... |
| Endotoxins | Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test. |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4. |
| Storage | Lyophilized protein should be stored at ≤ - 20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months. |
| Tag | C-Fc |
| Note | For research use only |
| Uniprot ID | AOA1BOGS01 |
| Expiration Date | 6 months from date of receipt. |