


Product Datasheet

Rel (Phospho-Ser503) Antibody (orb14754)

Description	Rel (Phospho-Ser503) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q04864
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of methanol-...