

Product Datasheet

Mouse BAFF Protein (orb1475402)

Description	Recombinant Mouse TNF Ligand Superfamily Member 13B is produced by our Mammalian expression system a...
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-mFc
Note	For research use only
Uniprot ID	Q9WU72
Expiration Date	6 months from date of receipt.