

Product Datasheet

Mouse CD27L Protein (orb1475422)

Description	Recombinant Mouse CD27 Ligand/CD70 is produced by our Mammalian expression system and the target gen...
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-10His
Note	For research use only
Uniprot ID	Q05A52
Expiration Date	6 months from date of receipt.