

Product Datasheet

Mouse CD7 Protein (orb1475428)

Description	Recombinant Mouse T-Cell Antigen CD7 is produced by our Mammalian expression system and the target g...
Endotoxins	Less than 0.01 EU/μg as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	P50283
Expiration Date	6 months from date of receipt.