

Product Datasheet

Human Siglec-8 Protein (orb1475431)

| | |
|------------------------|---|
| Description | Recombinant Human Sialic Acid-binding Ig-like Lectin 8 is produced by our Mammalian expression syste... |
| Endotoxins | Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test. |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of PBS,pH7.4. |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | C-mFc |
| Note | For research use only |
| Uniprot ID | Q9NYZ4 |
| Expiration Date | 6 months from date of receipt. |