

Product Datasheet

Human TSLP Protein (orb1475434)

Description	Recombinant Human Thymic stromal lymphopoietin is produced by our Mammalian expression system and th...
Endotoxins	Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-10His
Note	For research use only
Uniprot ID	Q969D9
Expiration Date	6 months from date of receipt.