



Product Datasheet

c-Jun (Phospho-Thr91) Antibody (orb14755)


Description	c-Jun (Phospho-Thr91) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05412
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of extracts from 2...



Western blot analysis of extract from He...