

Product Datasheet

Human NEDD8-conjugating enzyme UBE2F Protein (orb1475582)

| | |
|------------------------|--|
| Description | MADE TO ORDER Recombinant Human NEDD8-conjugating enzyme UBE2F is produced by E.coli expression sy... |
| Form/Appearance | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4 |
| Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Tag | C-His |
| Note | For research use only |
| Uniprot ID | Q969M7 |
| Expiration Date | 6 months from date of receipt. |