

## Product Datasheet

# Human Gamma-glutamylcyclotransferase Protein (orb1475584)

<b>Description</b>	MADE TO ORDER   Recombinant Human Gamma-glutamylcyclotransferase/GGCT is produced by E.coli expressi...
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
<b>Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>075223</b>
<b>Expiration Date</b>	6 months from date of receipt.