

## Product Datasheet

# Human Clusterin-associated protein 1 Protein (orb1475629)

<b>Description</b>	MADE TO ORDER   Recombinant Human Clusterin-associated protein 1/CLUAP1 is produced by E.coli expression system and the target gene encoding S2-F413 is expressed with a His tag at the C-terminus.
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
<b>Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q96AJ1</b>
<b>Expiration Date</b>	6 months from date of receipt.

---

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)