

Product Datasheet

Human Target of EGR1 protein 1 Protein (orb1475639)

Description	MADE TO ORDER Recombinant Human Target of EGR1 protein 1/TOE1 is produced by E.coli expression sys...
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze- thaw cycles.
Tag	C-His
Note	For research use only
Uniprot ID	Q96GM8
Expiration Date	6 months from date of receipt.