

Product Datasheet

Mouse IL-11 Protein (orb1475995)

Description	Recombinant Mouse IL-11 (N-Fc) is produced by our Mammalian expression system and the target gene en...
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-Fc-FLAG
Note	For research use only
Uniprot ID	P47873
Expiration Date	6 months from date of receipt.