


Product Datasheet


JunB (Phospho-Ser79) Antibody (orb14760)

Description	JunB (Phospho-Ser79) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P17275
Expiration Date	12 months from date of receipt.




P-Peptide

Immunohistochemical analysis of paraffin...



P-Peptide +

Immunohistochemical analysis of paraffin...



Immunofluorescence staining of methanol-
...