

Product Datasheet

Human Ubiquitin-conjugating enzyme E2 E2 Protein (orb1476212)

Description	MADE TO ORDER Recombinant Human Ubiquitin-conjugating enzyme E2 E2/UBE2E2 is produced by E.coli ex...
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Uniprot ID	Q96LR5
Expiration Date	6 months from date of receipt.