

Product Datasheet

Human AP-2 complex subunit alpha-2 Protein (orb1476329)

Description	MADE TO ORDER Recombinant Human AP-2 complex subunit alpha-2/AP2A2 is produced by E.coli expressio...
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Uniprot ID	O94973
Expiration Date	6 months from date of receipt.