

Product Datasheet

Human NEDD8 ultimate buster 1 Protein (orb1476588)

Description	MADE TO ORDER Recombinant Human NEDD8 ultimate buster 1/NUB1 is produced by E.coli expression syst...
Form/Appearance	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze- thaw cycles.
Tag	C-His
Note	For research use only
Uniprot ID	Q9Y5A7
Expiration Date	6 months from date of receipt.