

Product Datasheet

Mouse IL-2 Protein (orb1476620)

Description	Recombinant Mouse IL-2 is produced by our E.coli expression system and the target gene encoding Ala2...
Endotoxins	Less than 1 EU/ug
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Sodium Citrate, 0.2%Tween 80, pH 3.0
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Uniprot ID	P04351
Expiration Date	6 months from date of receipt.