


Product Datasheet


ATF2 (Phospho-Ser112 or 94) Antibody (orb14767)

Description	ATF2 (Phospho-Ser112 or 94) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P15336
Expiration Date	12 months from date of receipt.




P-Peptide

Immunohistochemical analysis of paraffin...



P-Peptide

Immunohistochemical analysis of paraffin...



Western blot analysis of ATF2(Phospho-Se...