

## Product Datasheet

### AHA1 antibody (orb147672)

**Description**

Rat monoclonal to Aha1 (HRP). Aha1 is a

is thought to st...

**Species/Host**

Rat

**Reactivity**

Human, Mouse, Rat

**Conjugation**

HRP

**Tested Applications**

ELISA, ICC, IHC, IP, WB

**Immunogen**

Recombinant Full Length Mouse Aha1 Protein

**Target**

AHA1

**Preservatives**

0.02M potassium phosphate, 0.15M sodium chloride, pH 7.2 \*Storage buffer may change when conjugated

**Concentration**

1.2 mg/ml

**Storage**

Conjugated antibodies should be stored at 4°C

**Note**

For research use only

**Application notes**

1 µg/ml of SMC-172 was sufficient for detection of Aha1 in 10 µg of rat tissue lysate by colorimetric immunoblot analysis using Goat anti-rat IgG:HRP as the secondary antibody.

**Isotype**

IgG2a Kappa

**Clonality**

Monoclonal

**MW**

38kDa

**Uniprot ID**

[Q8BK64](#)

**NCBI**

[NP\\_666148.1](#)

**Entrez**

[217737](#)

**Dilution Range**

WB (1:1000), IHC (1:100), ICC/IF (1:1000), IP (1:1000)

**Expiration Date**

12 months from date of receipt.