


Product Datasheet


Myc (Phospho-Thr358) Antibody (orb14769)

Description	Myc (Phospho-Thr358) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P01106
Expiration Date	12 months from date of receipt.




Immunofluorescence staining of methanol-

...



Immunohistochemical analysis of paraffin...



Western blot analysis of extracts from H...