



Product Datasheet

Elk1 (Phospho-Ser389) Antibody (orb14771)


| | |
|----------------------------|---|
| Description | Elk1 (Phospho-Ser389) Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P19419 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...



Western blot analysis of Elk1(Phospho-Se...