


Product Datasheet


GATA1 (Phospho-Ser142) Antibody (orb14774)

Description	GATA1 (Phospho-Ser142) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P15976
Expiration Date	12 months from date of receipt.




P-Peptide

Immunohistochemical analysis of paraffin...



P-Peptide

Immunohistochemical analysis of paraffin...



Western blot analysis of extracts from K...