


Product Datasheet


STAT4 (Phospho-Tyr693) Antibody (orb14779)

Description	STAT4 (Phospho-Tyr693) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q14765
Expiration Date	12 months from date of receipt.




P-Peptide -

Immunohistochemical analysis of paraffin...



P-Peptide +

Immunohistochemical analysis of paraffin...



Western blot analysis of extracts from H...