



Product Datasheet

STAT5a (phospho-Tyr694) antibody (orb14780)

Description	STAT5a (phospho-Tyr694) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human STAT5a around the phosphorylation site of Tyr694 (DG(p-Y)VK)
Target	phospho-STAT5a(Tyr694)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	91 kDa
Uniprot ID	P42229
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



IHC-P of rat lung tissue using STAT5a (p...