

---

## Product Datasheet

### Human BECN1 Protein (orb1478063)

**Description**

Recombinant Human Beclin-1(BECN1)

**Endotoxins**

Not test.

**Preservatives**

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

**Form/Appearance**

Liquid or Lyophilized powder

**Storage**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

**Tag**

N-terminal 10xHis-tagged

**Note**

For research use only

**Application notes**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Protein Sequence**

MEGSKTSNNSTMQVSFVCQRCSQPLKLDTSFKILDRVTIQELTAPLLTTAQAKPG  
 ETQEEETNSGEEPFIEIETPRQDGVSRRFIPPARMMSTESANSFTLIGEASDGGTME  
 NLSRRLKVTGDLFDIMSGQTDVDHPLCEECDTLLDQLDTQLNVTENECQNYKR  
 CLEILEQMNEEDDSEQLQMELELALAEERLIQELEDVEKNRKIVAENLEKVQAEAE  
 RLDQEEAQYQREYSEFKRQLELDDELKSVENQMRYAQTQLDKLKTNVFNAT  
 FHIWHSQGFGTINNFRGLRPSVPVEWNEINAAWGQTVLLLHALANKMGLKFQR  
 YRLVPYGNHSYLES�TDKSKELPLYCSGGLRFFWDNKFHAMVAFLDCVQQFKE  
 EVEKGETRFCLPYRMDVEKGIKIEDTGGSGGSYSIKTQFNSEEQWTKALKFMLTN  
 LKWGLAWVSSQFYNK

**Purity**

Greater than 85% as determined by SDS-PAGE.

**Source**

E.coli

**MW**

57.9 kDa

**Uniprot ID**
[Q14457](https://www.uniprot.org/entry/Q14457)
**Expiration Date**

6 months from date of receipt.