



Product Datasheet

Akt (Phospho-Ser473) Antibody (orb14785)

Description	Akt(Phospho-Ser473) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P31749
Expiration Date	12 months from date of receipt.




Western blot analysis of extract from He...



P-Peptide

Immunohistochemical analysis of paraffin...



P-Peptide +

IHC control with competing peptide