

Product Datasheet

Zika E (N2) antibody (orb1478821)

Description	Rabbit polyclonal antibody to Zika E (N2)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal aa 140-180 region of Zika virus (strain Mr 766) envelope protein E.
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	54
Uniprot ID	Q32ZE1
Expiration Date	12 months from date of receipt.