

Product Datasheet

Zika Pr (C) antibody (orb1478822)

Description	Rabbit polyclonal antibody to Zika Pr (C)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of Zika virus (strain Mr 766) peptide Pr
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	11
Uniprot ID	Q32ZE1
Expiration Date	12 months from date of receipt.