

## Product Datasheet

### Zika M (N) antibody (orb1478823)

<b>Description</b>	Rabbit polyclonal antibody to Zika M (N)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	N-terminal region of Zika virus (strain Mr 766) membrane protein M
<b>Preservatives</b>	0.01% NaN <sub>3</sub>
<b>Storage</b>	-20 °C
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	affinity purification
<b>MW</b>	9
<b>Uniprot ID</b>	<a href="#">Q32ZE1</a>
<b>Expiration Date</b>	12 months from date of receipt.