

Product Datasheet

ABI4 (C2) antibody (orb1478869)

Description	Rabbit polyclonal antibody to ABI4 (C2)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of arabidopsis thaliana ABI4
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	36
Uniprot ID	AOMES8
Expiration Date	12 months from date of receipt.