

## Product Datasheet

### EIN2 (N1) antibody (orb1478888)

<b>Description</b>	Rabbit polyclonal antibody to EIN2 (N1)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Plant
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	N-terminal region of arabidopsis thaliana EIN2
<b>Preservatives</b>	0.01% NaN <sub>3</sub>
<b>Storage</b>	-20 °C
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	affinity purification
<b>MW</b>	141
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/Q9S814">Q9S814</a>
<b>Expiration Date</b>	12 months from date of receipt.