


Product Datasheet

FLS2 (N) antibody (orb1478890)

Description	Rabbit polyclonal antibody to FLS2 (N)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of arabidopsis thaliana FLS2
Preservatives	0.01% NaN ₃
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	129
Uniprot ID	Q9FL28
Expiration Date	12 months from date of receipt.



Western blot analysis of total protein e...