

Product Datasheet

DCL4 (N2) antibody (orb1478989)

Description	Rabbit polyclonal antibody to DCL4 (N2)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of arabidopsis thaliana DCL4 (At5g20320)
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	191
Uniprot ID	P84634
Expiration Date	12 months from date of receipt.