

Product Datasheet

Mouse D2D ELISA Kit (orb1479384)

Description	Mouse D-Dimer (D2D) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Sample Type: plasma; competitive inhibition: Experimental Method
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	14.5pg/mL
Expiration Date	Please enquire.