

Product Datasheet

Mouse ADAM28 ELISA Kit (orb1479404)

Description	Mouse A Disintegrin And Metalloprotease 28 (ADAM28) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Sample Type: serum, plasma, tissue homogenates and other biological fluids; sandwich
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9JLN6
Sensitivity	0.058ng/mL
Entrez	13522
Expiration Date	Please enquire.