



Product Datasheet

CD329 antibody () (orb1480521)

Description	Mouse monoclonal antibody to CD329
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IP, WB
Immunogen	recombinant CD329-Hu IgGfC fusion protein
Target	CD329
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9Y336
Entrez	27180
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-fille...



Flow cytometry surface staining pattern ...