



Product Datasheet

NKp80 antibody () (orb1480524)

Description	Mouse monoclonal antibody to NKp80
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	recombinant human NKp80 extracellular domain
Target	NKp80
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9NZS2
Entrez	51348
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human NKp80 positive cells...



Flow cytometry surface staining pattern ...