


Product Datasheet


VEGFR2 (Phospho-Tyr1214) Antibody (orb14810)

Description	VEGFR2 (Phospho-Tyr1214) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P35968
Expiration Date	12 months from date of receipt.




P-Peptide -

Immunohistochemical analysis of paraffin...



P-Peptide +

Immunohistochemical analysis of paraffin...



Western blot analysis of extracts from S...