


Product Datasheet

Chk1 (Phospho-Ser280) Antibody (orb14859)

Description	Chk1 (Phospho-Ser280) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O14757
Expiration Date	12 months from date of receipt.



Western blot analysis of Chk1 (Phospho-Ser280) in HepG2 cells.