


Product Datasheet

Chk2 (Phospho-Ser516) Antibody (orb14864)

| | |
|----------------------------|---|
| Description | Chk2 (Phospho-Ser516) Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | O96017 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts from H...