
Product Datasheet

SHANK3 antibody (orb148738)

Description

Mouse monoclonal to SHANK3 (APC). Shank

identified th...

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	APC
Tested Applications	AM, ICC, IHC, IP, WB
Immunogen	Synthetic peptide amino acids 840-857 of rat Shank3
Target	SHANK3
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-336 was sufficient for detection of Shank3 in 10 µg COS cell lysate transiently transfected with Shank3 by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2b
Clonality	Monoclonal
MW	190kDa
Uniprot ID	Q9JLU4
NCBI	NP_067708.1
Entrez	59312
Dilution Range	WB (1:1000), IHC (1:100), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.