
Product Datasheet

GABA-A Receptor Delta antibody (orb148893)

Description

Mouse monoclonal to GABA-A R Delta

of the superfamily of fas...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

FITC

Tested

AM, ICC, WB

Applications
Immunogen

Synthetic peptide amino acids 15-34 (N-terminus) of rat GABA-A-R-Delta

Target

GABA A Receptor

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

2 µg/ml of SMC-345 was sufficient for detection of Delta1 GABA-A receptor in 10 µg of rat brain lysate by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2a

Clonality

Monoclonal

MW

55kDa

Uniprot ID

[P18506](#)

NCBI

[NP_058985.1](#)

Entrez

[29689](#)

Dilution Range

WB (1:1000), IHC (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.