



Product Datasheet

ASK1 (Phospho-Ser966) Antibody (orb14890)

Description	ASK1 (Phospho-Ser966) Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q99683
Expiration Date	12 months from date of receipt.



Western blot analysis of ASK1(Phospho-Ser...



Immunohistochemical analysis of formalin...