


Product Datasheet

HSP90B (Phospho-Ser254) Antibody (orb14901)

Description	HSP90B (Phospho-Ser254) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08238
Expiration Date	12 months from date of receipt.



Western blot analysis of HSP90B(Phospho-
...
...