


Product Datasheet

Synaptotagmin 1/2 (Phospho-Thr202/199) Antibody (orb14911)

Description	Synaptotagmin 1/2 (Phospho-Thr202/199) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P21579
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...